

Bion Multi Calibrator & QC

| Component | Bion Cal J501 EXP : 2026-12-12 | Bio-Rad Normal LOT : 89751 EXP : 2026-12-12 | | Bio-Rad Normal LOT : 89752 EXP : 2027-01-31 | |
|-------------------|--------------------------------------|---|------|---|-----|
| BION ALB | 3.6 | 4.0 | 0.2 | 2.7 | 0.4 |
| BION ALP | 240 | 136 | 10 | 517 | 25 |
| BION AMY | 185 | 66 | 8 | 378 | 40 |
| BION BUN | 48 | 16.0 | 1.5 | 49.4 | 3.5 |
| BION CHOL | 171 | 251 | 12 | 110 | 6 |
| BION CREA | 3.73 | 2.31 | 0.2 | 5.56 | 0.6 |
| BION D-BIL | 2.1 | 0.31 | 0.1 | 1.71 | 0.2 |
| BION GGT | 102 | 66 | 6.0 | 188 | 9 |
| BION Glucose | 212 | 87 | 5.0 | 293 | 15 |
| BION AST (GOT) | 116 | 46 | 4.0 | 236 | 10 |
| BION ALT (GPT) | 111 | 35 | 3.0 | 116 | 8.0 |
| BION HDL | 57 | 86 | 5.0 | 25 | 3.0 |
| BION LDH | 566 | 385 | 17.5 | 905 | 40 |
| BION LDL | 62 | 106 | 7.0 | 58 | 4.0 |
| BION T-BIL | 4.3 | 1.14 | 0.12 | 4.42 | 0.2 |
| BION TG | 145 | 194 | 10 | 87 | 5.0 |
| BION TP | 5.7 | 6.3 | 0.4 | 4.4 | 0.2 |
| BION UA | 5.3 | 4.55 | 0.35 | 9.7 | 0.7 |
| BION CA | 10.5 | 9.0 | 0.8 | 11.9 | 1.0 |
| BION P | 5.45 | 3.30 | 0.3 | 6.80 | 0.8 |
| BION MG | 2.9 | 2.2 | 0.3 | 4.7 | 0.4 |
| BION CK | 394 | 147 | 11 | 465 | 20 |

(주) 메타메디

약 30분간 실온에서 방치한 후 천천히
전도 혼합한 후에 사용합니다.

주소 : 대전광역시 서구 계백로 860번길 8, 203호

E-mail : mtmd@metamedi.net

Tel : 041-751-3826

Fax : 041-754-3827

