

# Bion Multi Calibrator & QC

Component	LiNEAR Multi CAL LOT : 19291 EXP : 2027-10-30	LiNEAR Normal LOT : 19292 EXP : 2027-10-30		LiNEAR Abormal LOT : 19293 EXP : 2027-10-30		Bio-Rad Normal LOT : 89751 EXP : 2027-01-31		Bio-Rad Abormal LOT : 89752 EXP : 2027-01-31	
BION ALB	3.1	3.7	0.32	2.8	0.23	3.8	0.32	2.7	0.23
BION ALP	358	170	15.80	367	28.10	118	10.1	439	40.3
BION AMY	315	103.0	7.72	318.0	23.80	78	6.4	456.5	41.8
BION BUN	52.8	20.6	1.80	53.7	3.96	14.2	1.2	42.4	3.9
BION CHOL	263	156	12.30	264	17.40	222	18.8	94.5	8.3
BION CREA	3.8	1.4	0.14	3.6	0.4	2.2	0.21	5.75	0.5
BION D-BIL	1.71	1.05	0.12	1.5	0.2	0.4	0.11	1.6	0.22
BION GGT	180	51.0	3.8	164.0	12.30	59	5.3	176.0	16.2
BION Glucose	275	110	8.25	266	21.00	84	7.1	284	24.1
BION AST (GOT)	140	39	3.98	136	13.90	39	3.3	201	19.0
BION ALT (GPT)	132	39	3.92	141	14.10	32	2.7	103	9.9
BION HDL	95	57	6.10	125	12.80	78.3	7.4	21.3	3.5
BION LDH	771	396	29.70	804	60.30	348	31.6	820	75.4
BION LDL	150	92	9.80	142	14.60	131.4	11.2	63.7	5.5
BION T-BIL	5.46	1.70	0.19	5.36	0.57	1.40	0.13	5.00	0.48
BION TG	260	88	7.22	260	21.10	171	14.5	77	6.5
BION TP	5.0	6.3	0.62	4.8	0.49	6.2	0.55	4.1	0.37
BION UA	9.8	6.1	0.57	9.8	0.81	4.9	0.42	10.3	0.9
BION CA	12.5	8.1	0.72	11.4	0.91	9.2	0.88	12.4	1.1
BION P	7.6	4.8	0.36	6.97	0.59	3.8	0.32	7.90	0.7
BION CK						137	7	450	23

(주) 메타메디

약 30분간 실온에서 방치한 후 천천히  
전도 혼합한 후에 사용합니다.

주소 : 대전광역시 서구 계백로 860번길 8, 203호

E-mail : mtmd@metamedi.net

Tel : 041-751-3826

Fax : 041-754-3827

**META-MEDI**  
(주) 메타메디

# Bion Multi Calibrator & QC

Component	Multi Cal (CFAS) LOT : 10759350 EXP : 2025-04-30	Multi Cal (CFAS) LOT : 76855401 EXP : 2026-03-31	Bio-Rad Normal LOT : 89751 EXP : 2027-01-31		Bio-Rad Abnormal LOT : 89752 EXP : 2027-01-31	
BION ALB	3.4	3.5	3.8	0.4	2.7	0.3
BION ALP	212	212	118	11.2	439	41.7
BION AMY	221	225	78	7.4	457	43.4
BION BUN	45.4	44.9	14.2	1.3	42.4	4.0
BION CHOL	145	147	222	21.0	94.5	9.0
BION CREA	3.7	3.58	2.2	0.2	5.75	0.5
BION D-BIL	2.2	2.4	0.4	0.12	1.6	0.48
BION GGT	97	97	59	5.6	176	16.7
BION Glucose	191	191	84	8.0	284	26.9
BION AST (GOT)	98	103	39	3.7	201	19.0
BION ALT (GPT)	108	106	32	3.0	103	9.8
BION HDL	45.2	48.8	78.3	7.4	21.3	2.0
BION LDH	494	538	348	33.1	820	77.9
BION LDL	83	85.5	131.4	12.5	63.7	6.0
BION T-BIL	4.7	5.2	1.40	0.1	5.00	0.5
BION TG	128	124	171	16.2	77	7.3
BION TP	5.1	5.2	6.2	0.6	4.1	0.4
BION UA	5.5	5.3	4.9	0.5	10.3	1.0
BION CA	10.4	10.7	9.2	0.9	12.4	1.2
BION P	6.1	6.0	3.8	0.4	7.90	0.8
BION MG	2.8	2.8	2.3	0.2	4.70	0.4
BION CK	387	365	137	13.0	450	42.8

(주) 메타메디

약 30분간 실온에서 방치한 후 천천히  
전도 혼합한 후에 사용합니다.

주소 : 대전광역시 서구 계백로 860번길 8, 203호

E-mail : mtmd@metamedi.net

Tel : 041-751-3826

Fax : 041-754-3827

**META-MEDI**  
(주) 메타메디